


AMINBIC
Advanced Materials Innovation



EXTRACTION KITS

Ronash Technology Pars

MagBic gDNA Extraction Kit

- **Cat. No.** D100100

- **Samples:** For blood samples like whole blood, plasma, serum and etc.

- **Apparatus and tools to be prepared:**

Magnetic Rack (MRCK100100).

- **Applications:**

PCR, RT-PCR, Quantitative Real-Time PCR, and Quantitative Real Time RT-PCR.



- **Description:**


MagBic gDNA Extraction Kit allows rapid and high yield extraction of ultra-pure DNA from various blood samples such as whole blood, plasma, serum and etc., using magnetic beads to efficiently minimize the isolation time.

- **Features and Benefits:**

- Minimized DNA extraction time from various samples using Magnetic Beads.
- No need of additional expensive instruments other than the kit and Magnetic Separation Rack.

- **Components:**

Components	Volume (ml)	Description
GDLB	50	Lysis Buffer
GDW	50	Washing Buffer
GDW I	100	Washing I Buffer
GDM	4.5	Magnetic Beads
GDE	10	Elution Buffer





Manufacturer: Ronash Technology Pars Co., Ltd.

Address: Unit 402, Third Floor, Electrochemical Institute, North Campus of Tehran University, Farshi
Moghadam st., North Kargar, Tehran, Iran

Tel: +989365650563

Post: 1449614535

Email: Info@aminbic.ir

Web: www.aminbic.ir